

Lab Settings

30/07/2025 4:21 pm BST

Summary

Lab Settings can be accessed via: Administration > Settings and selecting the 'Lab' category.

This page documents each setting.

Details

Configuration Settings

Categories

Account

Reporting

Diary

Reception

Treatment

Referral

Lab

Vet Scan Lab

Labelling

SMS

Email

DayLists

System

Miscellaneous

HTML Labelling

Smart Flow

VetCheck

Follow Up

GDPR

PACS

Message

Site

Esher SA

Save

Lab

External Lab

IDEXX Interlink Communication

Service Host

localhost

Service Port

4055

Client Port

4056

Client File Port

4057

File Download Timeout(seconds)

10

Connection Type

DIRECT_TCP

Server URL

http://localhost:8080/VoyagerEnterprise-Vc

Filtered Application Names (CSV)

For Idexx InterLink in Windows Terminal Session, please set INDIRECT TCP Connection Type and correct Server URL, example: http://[ServerNameOrIP][:PortOnlyIfRequired]/voyager-sessions-ejb/LabService?wsdl

HESKA Data Capture Utility(DCU)

DCU Output Folder Location :

C:\HeskaDcuResults

...

Please make sure the folder is shared (Read/Write Access) for all Merlin Users.

☒ Filter By Client Category

Idexx Interlink Communication.

These fields are used by MWI/Idexx to configure Idexx Interlink.

HESKA Data Capture Utility (DCU)

As above, these fields are used by MWI/Idexx to configure Idexx Interlink.

Filter by Client Category

Enabled = Allows users to filter Lab Results by Client Category in Utilities > General Lab.

Disabled = Do not allow users to filter Lab Results by Client Category in Utilities > General Lab.
