

Vet Scan Lab Settings


30/07/2025 4:37 pm BST

Summary

Vet Scan Settings can be accessed via: Administration > Settings and selecting the 'Vet Scan Lab' category.

This page is used to configure Abaxis Vet Scan Lab by our Engineers. This page documents each setting.

Details



Configuration Settings

Categories

- Account
- Reporting
- Diary
- Reception
- Treatment
- Referral
- Lab
- Vet Scan Lab**
- Labelling
- SMS
- Email
- DayLists
- System
- Miscellaneous
- HTML Labelling
- Smart Flow
- VetCheck
- Follow Up
- GDPR
- PACS
- Message

Site Esher SA

Vet Scan URL https://fuse03.abaxistest.com:8443/vetsync/v1

Vet Scan User pms

Vet Scan Password

Check results every 0 seconds

Vet Scan URL

The URL of the Vet Scan Fuse box.

Vet Scan User

Enter username here.

Vet Scan Password

Enter password here.

Check results every X seconds

Controls how often the system checks for results.