

Barcode Scanning

30/07/2025 1:00 pm BST

Summary


As an alternative to adding treatment using the methods documented in [Adding Treatment](#), items can be scanned using a hand-held device.

In order to use this functionality, [barcode data](#) will need to be imported into Merlin.


Details

Scan an Item - Add Treatment

1. Enter the 'Add Treatment' screen (Treatment > Add Treatment tab).
2. Select the 'Barcode' button.



3. The 'Scan Item' screen will display. Use your hand-held device to scan a product.

A screenshot of the 'Scan Item' screen. At the top, there are two buttons: 'Scanner Available' (highlighted with a dark blue background) and 'Scanner Not Available'. Below these is a section titled 'Details'. Inside the 'Details' section, there is a prompt 'Please scan item!' followed by a white rectangular input field. Below the input field, there is a message 'After scanning press Return to continue.' and a 'Cancel' button.

4. The item will add to the treatment draft.
5. Select 'Save' to save the item to treatment history.